



ATXN1L Antibody, Biotin conjugated

Product Code	CSB-PA002440LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0C7T5
Immunogen	Recombinant Human Ataxin-1-like protein (183-463AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ataxin-1-like (Brother of ataxin-1) (Brother of ATXN1), ATXN1L, BOAT BOAT1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ATXN1L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.