



BTK Antibody, FITC conjugated

Product Code	CSB-PA002867LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q06187
Immunogen	Recombinant Human Tyrosine-protein kinase BTK protein (254-443AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tyrosine-protein kinase BTK (EC 2.7.10.2) (Agammaglobulinemia tyrosine kinase) (ATK) (B-cell progenitor kinase) (BPK) (Bruton tyrosine kinase), BTK, AGMX1 ATK BPK
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	BTK
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.