




CCR5 Antibody

Product Code	CSB-PA004844LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51681
Immunogen	Recombinant Human C-C chemokine receptor type 5 protein (302-352AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, IF; Recommended dilution: IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	C-C chemokine receptor type 5 (C-C CKR-5) (CC-CKR-5) (CCR-5) (CCR5) (CHEMR13) (HIV-1 fusion coreceptor) (CD antigen CD195), CCR5, CMKBR5
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	CCR5

Image



Immunofluorescent analysis of HepG2 cells using CSB-PA004844LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)