



CD151 Antibody, HRP conjugated

Product Code	CSB-PA004880LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48509
Immunogen	Recombinant Human CD151 antigen protein (125-199AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	CD151 antigen (GP27) (Membrane glycoprotein SFA-1) (Platelet-endothelial tetraspan antigen 3) (PETA-3) (Tetraspanin-24) (Tspan-24) (CD antigen CD151), CD151, TSPAN24
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CD151
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.