



CD96 Antibody, FITC conjugated

Product Code	CSB-PA004971LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P40200
Immunogen	Recombinant Human T-cell surface protein tactile protein (306-395AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	T-cell surface protein tactile (Cell surface antigen CD96) (T cell-activated increased late expression protein) (CD antigen CD96), CD96
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CD96
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.