



CDH2 Antibody, HRP conjugated

Product Code	CSB-PA005045LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19022
Immunogen	Recombinant Human Cadherin-2 protein (746-906AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cadherin-2 (CDw325) (Neural cadherin) (N-cadherin) (CD antigen CD325), CDH2, CDHN NCAD
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CDH2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.