



CDH2 Antibody, Biotin conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA005045LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P19022 |
| Immunogen | Recombinant Human Cadherin-2 protein (746-906AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Cadherin-2 (CDw325) (Neural cadherin) (N-cadherin) (CD antigen CD325), CDH2, CDHN NCAD |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Signal Transduction |
| Target Names | CDH2 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |