




CHGB Antibody

Product Code	CSB-PA005345LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05060
Immunogen	Recombinant Human Secretogranin-1 protein (328-677AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Secretogranin-1 (Chromogranin-B) (CgB) (Secretogranin I) (Sgl) [Cleaved into: PE-11; GAWK peptide; CCB peptide], CHGB, SCG1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CHGB

Image



Western Blot

Positive WB detected in: HeLa whole cell lysate, HEK293 whole cell lysate, NIH/3T3 whole cell lysate

All lanes: CHGB antibody at 2.8µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 79 kDa

Observed band size: 79 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.