



CNR2 Antibody, HRP conjugated

Product Code	CSB-PA005679LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P34972
Immunogen	Recombinant Human Cannabinoid receptor 2 protein (305-360AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cannabinoid receptor 2 (CB-2) (CB2) (hCB2) (CX5), CNR2, CB2A CB2B
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	CNR2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.