



CORO1A Antibody, FITC conjugated

Product Code	CSB-PA005813LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31146
Immunogen	Recombinant Human Coronin-1A protein (363-432AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Coronin-1A (Coronin-like protein A) (Clipin-A) (Coronin-like protein p57) (Tryptophan aspartate-containing coat protein) (TACO), CORO1A, CORO1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CORO1A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.