




DUSP3 Antibody

Product Code	CSB-PA007255LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51452
Immunogen	Recombinant Human Dual specificity protein phosphatase 3 protein (14-96AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dual specificity protein phosphatase 3 (EC 3.1.3.16) (EC 3.1.3.48) (Dual specificity protein phosphatase VHR) (Vaccinia H1-related phosphatase) (VHR), DUSP3, VHR
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DUSP3

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA007255LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.