



F9 Antibody, Biotin conjugated

Product Code	CSB-PA007936LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00740
Immunogen	Recombinant Human Coagulation factor IX protein (144-239AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Coagulation factor IX (EC 3.4.21.22) (Christmas factor) (Plasma thromboplastin component) (PTC) [Cleaved into: Coagulation factor IXa light chain; Coagulation factor IXa heavy chain], F9
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	F9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.