




FGFR1 Antibody

Product Code	CSB-PA008642LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11362
Immunogen	Recombinant Human Fibroblast growth factor receptor 1 protein (22-96AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IF; Recommended dilution: WB:1:2000-1:10000, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fibroblast growth factor receptor 1 (FGFR-1) (EC 2.7.10.1) (Basic fibroblast growth factor receptor 1) (BFGFR) (bFGF-R-1) (Fms-like tyrosine kinase 2) (FLT-2) (N-sam) (Proto-oncogene c-Fgr) (CD antigen CD331), FGFR1, BFGFR CEK FGFBR FLG FLT2 HBGFR
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	FGFR1

Image



Western Blot

Positive WB detected in: A549 whole cell lysate

All lanes: FGFR1 antibody at 3.4µg/ml


Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 92, 74, 7, 83, 64, 82, 75,

56, 17, 34, 91, 96 kDa

Observed band size: 92 kDa



Immunofluorescent analysis of HeLa cells using CSB-PA008642LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.