



FGFR1 Antibody, FITC conjugated

Product Code	CSB-PA008642LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11362
Immunogen	Recombinant Human Fibroblast growth factor receptor 1 protein (22-96AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fibroblast growth factor receptor 1 (FGFR-1) (EC 2.7.10.1) (Basic fibroblast growth factor receptor 1) (BFGFR) (bFGF-R-1) (Fms-like tyrosine kinase 2) (FLT-2) (N-sam) (Proto-oncogene c-Fgr) (CD antigen CD331), FGFR1, BFGFR CEK FGFB R FLG FLT2 HBGFR
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	FGFR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.