



GRM1 Antibody, Biotin conjugated

Product Code	CSB-PA009931LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13255
Immunogen	Recombinant Human Metabotropic glutamate receptor 1 protein (880-1017AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Metabotropic glutamate receptor 1 (mGluR1), GRM1, GPRC1A MGLUR1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GRM1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.