



PSMB5 Antibody, HRP conjugated

Product Code	CSB-PA010083LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P28074
Immunogen	Recombinant Human Proteasome subunit beta type-5 protein (106-188AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Proteasome subunit beta type-5 (EC 3.4.25.1) (Macropain epsilon chain) (Multicatalytic endopeptidase complex epsilon chain) (Proteasome chain 6) (Proteasome epsilon chain) (Proteasome subunit MB1) (Proteasome subunit X), PSMB5, LMPX MB1 X
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PSMB5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.