



HNF4A Antibody, Biotin conjugated

Product Code	CSB-PA010594LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41235
Immunogen	Recombinant Human Hepatocyte nuclear factor 4-alpha protein (130-330AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hepatocyte nuclear factor 4-alpha (HNF-4-alpha) (Nuclear receptor subfamily 2 group A member 1) (Transcription factor 14) (TCF-14) (Transcription factor HNF-4), HNF4A, HNF4 NR2A1 TCF14
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	HNF4A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.