



HPSE Antibody, Biotin conjugated

Product Code	CSB-PA010716LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y251
Immunogen	Recombinant Human Heparanase protein (263-523AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Heparanase (EC 3.2.1.166) (Endo-glucuronidase) (Heparanase-1) (Hpa1) [Cleaved into: Heparanase 8 kDa subunit; Heparanase 50 kDa subunit], HPSE, HEP HPA HPA1 HPR1 HPSE1 HSE1
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	HPSE
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.