



IER3 Antibody, FITC conjugated

Product Code	CSB-PA011002LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P46695
Immunogen	Recombinant Human Radiation-inducible immediate-early gene IEX-1 protein (1-82AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Radiation-inducible immediate-early gene IEX-1 (Differentiation-dependent gene 2 protein) (Protein DIF-2) (Immediate early protein GLY96) (Immediate early response 3 protein) (PACAP-responsive gene 1 protein) (Protein PRG1), IER3, DIF2 IEX1 PRG1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	IER3