



IER3 Antibody, FITC conjugated

Product Code	CSB-PA011002LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P46695
Immunogen	Recombinant Human Radiation-inducible immediate-early gene IEX-1 protein (1-82AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Radiation-inducible immediate-early gene IEX-1 (Differentiation-dependent gene 2 protein) (Protein DIF-2) (Immediate early protein GLY96) (Immediate early response 3 protein) (PACAP-responsive gene 1 protein) (Protein PRG1), IER3, DIF2 IEX1 PRG1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	IER3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.