



IGF2BP3 Antibody, Biotin conjugated

Product Code	CSB-PA011092LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00425
Immunogen	Recombinant Human Insulin-like growth factor 2 mRNA-binding protein 3 protein (1-100AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Insulin-like growth factor 2 mRNA-binding protein 3 (IGF2 mRNA-binding protein 3) (IMP-3) (IGF-II mRNA-binding protein 3) (KH domain-containing protein overexpressed in cancer) (hKOC) (VICKZ family member 3), IGF2BP3, IMP3 KOC1 VICKZ3
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	IGF2BP3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.