



IL4 Antibody, Biotin conjugated

Product Code	CSB-PA011659YD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05112
Immunogen	Recombinant Human Interleukin-4 protein (25-153AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-4 (IL-4) (B-cell stimulatory factor 1) (BSF-1) (Binetrakin) (Lymphocyte stimulatory factor 1) (Pitrakinra), IL4
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IL4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.