



# KCNE1 Antibody, HRP conjugated

Product Code	CSB-PA012025LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15382
Immunogen	Recombinant Human Potassium voltage-gated channel subfamily E member 1 protein (85-129AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Potassium voltage-gated channel subfamily E member 1 (Delayed rectifier potassium channel subunit Isk) (IKs producing slow voltage-gated potassium channel subunit beta Mink) (Minimal potassium channel), KCNE1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	KCNE1