



LANCL1 Antibody, HRP conjugated

Product Code	CSB-PA012742LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43813
Immunogen	Recombinant Human LanC-like protein 1 protein (50-150AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	LanC-like protein 1 (40 kDa erythrocyte membrane protein) (p40), LANCL1, GPR69A
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	LANCL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.