



LIPE Antibody, Biotin conjugated

Product Code	CSB-PA012975LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q05469
Immunogen	Recombinant Human Hormone-sensitive lipase protein (1-251AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hormone-sensitive lipase (HSL) (EC 3.1.1.79), LIPE
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	LIPE
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.