



MAP2K4 Antibody, Biotin conjugated

Product Code	CSB-PA013413LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P45985
Immunogen	Recombinant Human Dual specificity mitogen-activated protein kinase kinase 4 protein (259-398AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dual specificity mitogen-activated protein kinase kinase 4 (MAP kinase kinase 4) (MAPKK 4) (EC 2.7.12.2) (JNK-activating kinase 1) (MAPK/ERK kinase 4) (MEK 4) (SAPK/ERK kinase 1) (SEK1) (Stress-activated protein kinase kinase 1) (SAPK kinase 1) (SAPKK-1) (SAPKK1) (c-Jun N-terminal kinase kinase 1) (JNKK), MAP2K4, JNKK1 MEK4 MKK4 PRKMK4 SEK1 SERK1 SKK1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MAP2K4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.