




MDM4 Antibody

Product Code	CSB-PA013627LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15151
Immunogen	Recombinant Human Protein Mdm4 protein (208-306AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein Mdm4 (Double minute 4 protein) (Mdm2-like p53-binding protein) (Protein Mdmx) (p53-binding protein Mdm4), MDM4, MDMX
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MDM4

Image



Western Blot

Positive WB detected in: HeLa whole cell lysate, HepG2 whole cell lysate, Mouse stomach tissue

All lanes: MDM4 antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 55, 15, 14, 19, 50 kDa

Observed band size: 55 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.