



MDM4 Antibody, FITC conjugated

Product Code	CSB-PA013627LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15151
Immunogen	Recombinant Human Protein Mdm4 protein (208-306AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein Mdm4 (Double minute 4 protein) (Mdm2-like p53-binding protein) (Protein Mdmx) (p53-binding protein Mdm4), MDM4, MDMX
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MDM4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.