



MFN2 Antibody, Biotin conjugated

Product Code	CSB-PA013756LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95140
Immunogen	Recombinant Human Mitofusin-2 protein (61-158AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitofusin-2 (EC 3.6.5.-) (Transmembrane GTPase MFN2), MFN2, CPRP1 KIAA0214
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	MFN2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.