



MX1 Antibody, Biotin conjugated

Product Code	CSB-PA015249LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20591
Immunogen	Recombinant Human Interferon-induced GTP-binding protein Mx1 protein (7-172AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interferon-induced GTP-binding protein Mx1 (Interferon-induced protein p78) (IFI-78K) (Interferon-regulated resistance GTP-binding protein MxA) (Myxoma resistance protein 1) (Myxovirus resistance protein 1) [Cleaved into: Interferon-induced GTP-binding protein Mx1, N-terminally processed], MX1
Immunogen Species	Homo sapiens (Human)
Research Area	Microbiology
Target Names	MX1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.