





MYBL2 Antibody

Product Code	CSB-PA015264LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10244
Immunogen	Recombinant Human Myb-related protein B protein (188-366AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myb-related protein B (B-Myb) (Myb-like protein 2), MYBL2, BMYB
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MYBL2

Image



Immunohistochemistry of paraffin-embedded human pancreatic cancer using CSB-PA015264LA01HU at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA015264LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)