



NNMT Antibody, FITC conjugated

Product Code	CSB-PA015906LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P40261
Immunogen	Recombinant Human Nicotinamide N-methyltransferase protein (86-157AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nicotinamide N-methyltransferase (EC 2.1.1.1), NNMT
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NNMT