




PAFAH1B1 Antibody

Product Code	CSB-PA017383LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43034
Immunogen	Recombinant Human Platelet-activating factor acetylhydrolase IB subunit alpha protein (40-155AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IF; Recommended dilution: IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Platelet-activating factor acetylhydrolase IB subunit alpha (Lissencephaly-1 protein) (LIS-1) (PAF acetylhydrolase 45 kDa subunit) (PAF-AH 45 kDa subunit) (PAF-AH alpha) (PAFAH alpha), PAFAH1B1, LIS1 MDCR MDS PAFAHA
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PAFAH1B1

Image



Immunofluorescent analysis of HeLa cells using CSB-PA017383LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.