



PAFAH1B1 Antibody, HRP conjugated

Product Code	CSB-PA017383LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43034
Immunogen	Recombinant Human Platelet-activating factor acetylhydrolase IB subunit alpha protein (40-155AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Platelet-activating factor acetylhydrolase IB subunit alpha (Lissencephaly-1 protein) (LIS-1) (PAF acetylhydrolase 45 kDa subunit) (PAF-AH 45 kDa subunit) (PAF-AH alpha) (PAFAH alpha), PAFAH1B1, LIS1 MDCR MDS PAFAHA
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PAFAH1B1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.