



# PAK3 Antibody, Biotin conjugated

Product Code	CSB-PA017407LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75914
Immunogen	Recombinant Human Serine/threonine-protein kinase PAK 3 protein (70-190AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase PAK 3 (EC 2.7.11.1) (Beta-PAK) (Oligophrenin-3) (p21-activated kinase 3) (PAK-3), PAK3, OPHN3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PAK3