




PER1 Antibody

Product Code	CSB-PA017786LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15534
Immunogen	Recombinant Human Period circadian protein homolog 1 protein (1-294AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Period circadian protein homolog 1 (hPER1) (Circadian clock protein PERIOD 1) (Circadian pacemaker protein Rigi), PER1, KIAA0482 PER RIGUI
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PER1

Image



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using CSB-PA017786LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.