



PER1 Antibody, FITC conjugated

Product Code	CSB-PA017786LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15534
Immunogen	Recombinant Human Period circadian protein homolog 1 protein (1-294AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Period circadian protein homolog 1 (hPER1) (Circadian clock protein PERIOD 1) (Circadian pacemaker protein Rigiui), PER1, KIAA0482 PER RIGUI
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PER1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.