



PRSS2 Antibody, Biotin conjugated

Product Code	CSB-PA018814LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07478
Immunogen	Recombinant Human Trypsin-2 protein (24-247AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Trypsin-2 (EC 3.4.21.4) (Anionic trypsinogen) (Serine protease 2) (Trypsin II), PRSS2, TRY2 TRYP2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PRSS2