



PTGR1 Antibody, FITC conjugated

Product Code	CSB-PA018983LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14914
Immunogen	Recombinant Human Prostaglandin reductase 1 protein (202-282AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Prostaglandin reductase 1 (PRG-1) (EC 1.3.1.-) (15-oxoprostaglandin 13-reductase) (EC 1.3.1.48) (NADP-dependent leukotriene B4 12-hydroxydehydrogenase) (EC 1.3.1.74), PTGR1, LTB4DH
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	PTGR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.