



PTK2 Antibody, HRP conjugated

Product Code	CSB-PA018994LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q05397
Immunogen	Recombinant Human Focal adhesion kinase 1 protein (363-464AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Focal adhesion kinase 1 (FADK 1) (EC 2.7.10.2) (Focal adhesion kinase-related nonkinase) (FRNK) (Protein phosphatase 1 regulatory subunit 71) (PPP1R71) (Protein-tyrosine kinase 2) (p125FAK) (pp125FAK), PTK2, FAK FAK1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	PTK2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.