



RPS6KA4 Antibody, Biotin conjugated

Product Code	CSB-PA020467LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75676
Immunogen	Recombinant Human Ribosomal protein S6 kinase alpha-4 protein (237-393AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ribosomal protein S6 kinase alpha-4 (S6K-alpha-4) (EC 2.7.11.1) (90 kDa ribosomal protein S6 kinase 4) (Nuclear mitogen- and stress-activated protein kinase 2) (Ribosomal protein kinase B) (RSKB), RPS6KA4, MSK2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RPS6KA4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.