



SFRP1 Antibody, HRP conjugated

Product Code	CSB-PA021138LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N474
Immunogen	Recombinant Human Secreted frizzled-related protein 1 protein (170-240AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Secreted frizzled-related protein 1 (FRP-1) (sFRP-1) (Secreted apoptosis-related protein 2) (SARP-2), SFRP1, FRP FRP1 SARP2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SFRP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.