



SFRP4 Antibody, Biotin conjugated

Product Code	CSB-PA021140LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6FHJ7
Immunogen	Recombinant Human Secreted frizzled-related protein 4 protein (293-346AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Secreted frizzled-related protein 4 (sFRP-4) (Frizzled protein, human endometrium) (FrpHE), SFRP4, FRPHE
Immunogen Species	Homo sapiens (Human)
Research Area	Stem Cells
Target Names	SFRP4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.