



SIK1 Antibody, HRP conjugated

Product Code	CSB-PA021321LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P57059
Immunogen	Recombinant Human Serine/threonine-protein kinase SIK1 protein (492-750AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase SIK1 (EC 2.7.11.1) (Salt-inducible kinase 1) (SIK-1) (Serine/threonine-protein kinase SNF1-like kinase 1) (Serine/threonine-protein kinase SNF1LK), SIK1, SIK SNF1LK
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SIK1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.