



SIK1 Antibody, FITC conjugated

Product Code	CSB-PA021321LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P57059
Immunogen	Recombinant Human Serine/threonine-protein kinase SIK1 protein (492-750AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase SIK1 (EC 2.7.11.1) (Salt-inducible kinase 1) (SIK-1) (Serine/threonine-protein kinase SNF1-like kinase 1) (Serine/threonine-protein kinase SNF1LK), SIK1, SIK SNF1LK
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SIK1