



SLC6A3 Antibody, FITC conjugated

Product Code	CSB-PA021701LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01959
Immunogen	Recombinant Human Sodium-dependent dopamine transporter protein (161-237AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sodium-dependent dopamine transporter (DA transporter) (DAT) (Solute carrier family 6 member 3), SLC6A3, DAT1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SLC6A3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.