



STIM2 Antibody, Biotin conjugated

Product Code	CSB-PA022830LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P246
Immunogen	Recombinant Human Stromal interaction molecule 2 protein (471-656AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Stromal interaction molecule 2, STIM2, KIAA1482
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	STIM2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.