



TANK Antibody, Biotin conjugated

Product Code	CSB-PA023116LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92844
Immunogen	Recombinant Human TRAF family member-associated NF-kappa-B activator protein (4-98AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	TRAF family member-associated NF-kappa-B activator (TRAF-interacting protein) (I-TRAF), TANK, ITRAF TRAF2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TANK
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.