



Tgfb1 Antibody, HRP conjugated

Product Code	CSB-PA023446LB01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04202
Immunogen	Recombinant Mouse Transforming growth factor beta-1 protein (30-278AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transforming growth factor beta-1 (TGF-beta-1) [Cleaved into: Latency-associated peptide (LAP)], Tgfb1
Immunogen Species	Mus musculus (Mouse)
Research Area	Cardiovascular
Target Names	Tgfb1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.