



Tgfb1 Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA023446LC01MO |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P04202 |
| Immunogen | Recombinant Mouse Transforming growth factor beta-1 protein (30-278AA) |
| Raised In | Rabbit |
| Species Reactivity | Mouse |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Transforming growth factor beta-1 (TGF-beta-1) [Cleaved into: Latency-associated peptide (LAP)], Tgfb1 |
| Immunogen Species | Mus musculus (Mouse) |
| Research Area | Cardiovascular |
| Target Names | Tgfb1 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |